
TRICVALVE®

TRANSCATHETER BICAVAL VALVE WITH DELIVERY SYSTEM

EN TricValve® Transcatheter Bicaval Valve with Delivery System | **DE** TricValve® Bikavale Transkatheterklappe mit Einführsystem | **FR** TricValve® Valve Bicavale Transcathéter avec Système de mise en place | **ES** TricValve® Válvula Transcatéter Bicaval con Sistema de Liberación | **IT** TricValve® Valvola Bicaval Transcatetere con Sistema di Rilascio | **PT** TricValve® Válvula Bicaval Transcateter com Sistema de Implantação

EN INSTRUCTIONS FOR USE | **DE** GEBRAUCHSANWEISUNG | **FR** MODE D'EMPLOI | **ES** INSTRUCCIONES DE USO | **IT** ISTRUZIONI PER L'USO | **PT** INSTRUÇÕES DE USO

CE 0124

P&F PRODUCTS & FEATURES

GmbH

**EN CONTENT | DE INHALT | FR CONTENU |
ES CONTENIDO | IT CONTENUTO | PT CONTEÚDO**

| | |
|--------------------------------------|----|
| EN INSTRUCTIONS FOR USE..... | 4 |
| DE GEBRAUCHSANWEISUNG | 20 |
| FR MODE D'EMPLOI | 34 |
| ES INSTRUCCIONES DE USO | 49 |
| IT ISTRUZIONI PER L'USO | 64 |
| PT INSTRUÇÕES DE USO | 79 |

TricValve® Transcatheter Bicaval Valve with Delivery System

CAUTION: Implantation of TricValve® Transcatheter Bicaval Valves should only be performed by physicians who have received or undergone P+F Products + Features GmbH training on usage of device. The devices are supplied sterile for single use only. After use, dispose the TricValve® Delivery System in accordance with local regulations and hospital procedures. Do not re-sterilize.

STERILE: The TricValve® Transcatheter Bicaval Valves System is supplied in sterilized condition with Ethylene oxide (EtO) gas.

1. PRODUCT DESCRIPTION

The TricValve® Transcatheter Bicaval Valves System is a Bicaval Transcatheter tricuspid valve implantation system, which includes the TricValve® Transcatheter Bicaval Valve for Superior Vena Cava (SVC), the TricValve® Transcatheter Bicaval Valve for Inferior Vena Cava (IVC). The TricValve® Transcatheter Bicaval Valves are already pre-mounted into the TricValve® Delivery System.

1.1. TRICVALVE® TRANSCATHETER BICAVAL VALVE (BIOPROSTHESIS)

The bioprosthesis is developed for treatment of caval reflux present in cases of severe tricuspid regurgitation with a high risk for open-heart surgery. The bioprosthesis is designed to treat severe tricuspid regurgitation without removal of the defective tricuspid valve. The bioprosthesis is available in two different diameters for each models (SVC and IVC) specifically designed to adapt to the anatomic features of the superior and inferior vena cava (Figure 1 and Table 1).

The bioprosthesis consists of a tubular metallic structure of nitinol which is self-expandable and radiopaque with three valve leaflets of bovine pericardium sutured and complemented by a skirt of polyester to avoid paravalvular leaks. The bioprosthesis leaflets are processed with anti-calcification treatment as well as chemical dehydration.

The valve is available in different diameters as shown in Table 1.


Figure 1:
TricValve® Transcatheter Bicaval Valves for Superior Vena Cava SVC (a) and Inferior Vena Cava IVC (b)

| TRICVALVE® MODEL | VALVE SIZE | PROXIMAL DIAMETER | DISTAL DIAMETER | LENGTH AFTER DEPLOYMENT |
|------------------|------------|-------------------|-----------------|-------------------------|
| SVC 25 | 25 | 25 | 20 | 66,60 |
| SVC 29 | 29 | 29 | 20 | 69,10 |
| IVC 31 | 31 | 34 | 38 | 65 |
| IVC 35 | 35 | 38 | 47 | 65 |

Table 1: Sizes of TricValve® Transcatheter Bicaval Valves (mm)


1.2. TRICVALVE® DELIVERY SYSTEM (CATHETER)

The catheter (Figure 2) is an enclosure retraction type system and used to deploy the bioprosthesis. It has anchors in its system for safe and accurate bioprosthesis deployment. The distal end of the system has an atraumatic radiopaque tip and protective sheath. A capsule at distal end covers and maintains the bioprosthesis in a crimped position. A stabilizer tube is fixed at the handle and extends outside the catheter shaft. It provides a barrier between the inner catheter shaft and vessel walls, thus enabling the catheter to retract freely. The delivery system is compatible with 0.889mm (0.035 inch) guide wire.

The handle includes a macro slider for opening and closing the bioprosthesis housing and a micro adjusting knob to facilitate accurate release of the bioprosthesis. The micro knob rotates clockwise to open the housing and in anticlockwise direction for close the capsule. The delivery system has a flush port which is used to hydrate the bioprosthesis leaflets and remove air before usage.

There is a single catheter model which is used for delivery of both bioprosthesis models.

Figure 2: TricValve® Delivery System

2. WARNINGS AND PRECAUTIONS

Before use, carefully read all warnings, precautions, and instructions for use for all components of the system. Failure to read and follow all instructions or failure to observe all stated warnings/instructions could cause serious injury or death to the patient.

2.1. WARNINGS

The devices are designed for single use only. Do not re-sterilize or reuse the devices. There is no data to support the sterility, non-pyrogenic, and functionality of the devices after reprocessing or re-sterilisation.

- Verify that relevant patient anatomical parameters are suitable and within the specifications for performing the procedure.
- Incorrect sizing of the bioprosthesis may lead to paravalvular leak, migration, embolization and/or cardio-vascular rupture.
- Accelerated deterioration of the bioprosthesis may occur in patients with an altered calcium metabolism.
- Valve leaflets mishandled or damaged during any part of the procedure will require replacement of the bioprosthesis.
- Caution should be exercised in implanting a valve in patients with clinically significant coronary artery disease.
- Caution should be exercised in patients with relevant left sided valve disease.
- Do not use the system if the tamper evident seal is broken or the valve is damaged, or the expiration date has elapsed.
- Do not use the catheter if the packaging or any components are not sterile, have been opened or are damaged (e.g. kinked or stretched).
- Follow the protocol for measuring relevant anatomical parameters of the patient before selecting the bioprosthesis for the treatment as described in the section "4. PATIENT SELECTION AND TREATMENT".
- After the delivery system was inserted into a patient, do not attempt to reload it in the same or another release device. Dispose the delivery system and catheter; do not attempt to reuse any component.
- Do not open the contents of the package until you are sure about the deployment and the appropriate size of the bioprosthesis.
- Do not handle or manipulate the bioprosthesis with sharp or pointed objects.
- Do not use the delivery system with any damage.
- Any mechanical failures of the delivery system may result in complications and severe damage to the patient.
- In order to avoid contamination of the delivery system, do not use gloves with powder.
- Use the products before the use by-date.
- Store the bioprosthesis between 8 and 25°C (46°F to 77°F).
- Carefully remove the catheter from the packaging to prevent damage to the catheter.

2.2. PRECAUTIONS

Do not use the TricValve® Transcatheter Bicaval Valves System

- If the sterile package is damaged;
- If the bioprosthesis is damaged;
- If the delivery system is damaged;
- If the delivery system is unable to flush;
- If the use-by-date has elapsed.

2.3. REPOSITIONING PRECAUTIONS

- After initiation of the release, withdrawal of the bioprosthesis (e.g. catheter use) is not recommended. A withdrawal of a partially released bioprosthesis can cause mechanical failure of the releasing device, damage to a vena cava, an iliac or femoral vein, as well as vascular complications, prosthetic valve dysfunction (including an abnormal position of the device), embolization, atrium and / or lead to emergency surgery.
- During release, the bioprosthesis may be retracted, if the vena cava has not been touched. After contact with the vein, it will not be possible to retract the bioprosthesis.
- After release, repositioning of the bioprosthesis (e.g. using tweezers or loops) is not recommended. Repositioning a released bioprosthesis may damage the vena cava or the atrium, as well as lead to vascular complications, prosthetic valve dysfunction (including an abnormal position of the device), embolization, and / or lead to emergency surgery.
- Do not try to remove the delivery system if some piece of the bioprosthesis is still connected with the delivery system. Ensure that there was a full bioprosthesis release prior to initiating the delivery system withdrawal procedure.

2.4. PRECAUTIONS AFTER IMPLANTATION

- Make sure the capsule is closed before removing the delivery system. If you feel an increase in resistance when removing the delivery system, do not force the passage. Increased resistance may be indicative of a problem and may damage the device and / or cause injury to the patient if the passage is forced. If the cause of the resistance cannot be determined or corrected, slowly rotate the delivery system and inspect whether the delivery system is complete or if there is any bent or ruptured region. One must always make sure that there is no connection of the delivery system with the valve during the removal of it.
- Care should be taken during removal of the delivery system from the patient, maintaining radioscopic visualization.
- Care should be taken when crossing a TricValve® Transcatheter Bicaval Valve implanted with ancillary devices.

3. POTENTIAL ADVERSE EVENTS

Potential risks associated with the implantation of the bioprosthesis may include, but are not limited to, the following:

- Abnormal laboratory values (including electrolyte imbalance);
- Allergic reaction to anticoagulation agents, contrast medium or anesthesia;
- Anemia;
- Angina;
- Asymptomatic cerebral infarct;
- Cardiac arrest;
- Cardiac arrhythmias;
- Cardiac insufficiency;
- Cardiogenic shock;
- Cardiovascular injury (including rupture, perforation or dissection of vessels, ventricles, atrium, or cardiac structures that may require intervention);
- Complications related to vascular access (eg, dissection, perforation, pain, hemorrhage, hematoma, pseudo aneurysm, irreversible nerve damage, compartment syndrome, arteriovenous fistula, stenosis);
- Death;
- Dyspnea;
- Embolization of the auxiliary device;
- Encephalopathy;
- Erosion of tissues;
- Fever;
- Heart failure or low cardiac output;
- Heart murmur;
- Hemolysis;

- Hypotension or hypertension;
- Infection (including septicemia and endocarditis);
- Inflammation;
- Injury of the vena cava;
- Intestinal ischemia;
- Migraine, transient ischemic attack (TIA) or other neurological deficits;
- Multiple organ insufficiency;
- Myocardial infarction;
- Myocardial ischemia;
- New non-emergency surgery;
- Pericardial effusion;
- Peripheral ischemia;
- Permanent disability;
- Pleural effusion;
- Prosthetic valve dysfunction including, but not limited to, fracture; bending (deformed conformation) of the bioprosthesis structure; sub expansion of the valve structure; calcification; pannus; wear, tear, prolapse or retraction in leaflets; poor valve coaptation; failure or disruption of suture; leakage; incorrect sizing (prosthesis-patient maladaptation); incorrect positioning (too high or too low) / incorrect placement; regurgitation and stenosis;
- Pulmonary edema;
- Renal insufficiency or renal failure (including acute renal injury);
- Respiratory insufficiency or respiratory failure;
- Severe or small bleeding that may or may not require transfusion or intervention (including fatal hemorrhage or disabling bleeding);
- Syncope;
- Thrombosis / clot (including valve thrombosis);
- Valve migration / valve embolization.

4. PATIENT SELECTION AND TREATMENT

4.1. INDICATIONS

TricValve® Transcatheter Bicaval Valves System is a system of two self-expanding biological valves for the treatment of patients with hemodynamically relevant tricuspid insufficiency and caval reflux. The bioprosthesis are implanted percutaneously into the inferior and superior vena cava without disturbing the native tricuspid valve. It is especially intended for use in patients at extreme risk for open heart surgery or who are deemed inoperable by a Heart Team.

4.2. CONTRAINDICATIONS

The TricValve® Transcatheter Bicaval Valves System is contraindicated in patients who have any of the following conditions:

- Evolutionary or recent stroke;
- Cerebrovascular accident (CVA) evolutionary or recent;
- Recent myocardial infarction (<30 days);
- Known hypersensitivity, allergy or contraindication to device's components, e.g. nitinol;
- Known hypersensitivity to vitamin K antagonists, heparin and other oral anti-coagulants, or sensitivity to contrast medium that cannot be adequately premedicated;
- Sepsis, including active endocarditis;
- Thrombosis of the lower venous system or vena cava filter;
- Known hypersensitivity to vitamin K antagonists, heparin and other oral anti-coagulants, or sensitivity to contrast medium that cannot be adequately premedicated;
- The patient must be 18 years of age or older;
- Creatinine clearance <20 ml / min;
- Vascular conditions (e.g. stenosis, tortuosity) that make insertion and endovascular access impossible to the upper and lower vena cava;
- Bleeding diathesis or coagulopathy or patient refusing blood transfusion;
- Active gastritis or peptic ulcer;
- Pregnancy.

4.3. TREATMENT CONSIDERATIONS

The risks and benefits previously described should be carefully considered for each patient before using TricValve® Transcatheter Bicaval Valves System. The risks of long-term anticoagulant and/or antiplatelet therapy should be taken into account. Generally, implantation of a bioprosthesis must be avoided in patients with any of the following conditions:

- High risk of bleeding (e.g. patients with active gastritis or recent peptic ulcer); and
- Pre-existing conditions of morbidity that increase the risk of poor initial results (e.g. diabetes mellitus, renal failure or severe obesity).

4.4. VALVE-IN-VALVE CONSIDERATIONS

The safety and efficacy of implantation of a second bioprosthesis within the initial bioprosthesis has not been studied. However, if a second bioprosthesis is needed within the initial bioprosthesis to improve valve function, valve size and patient anatomy should be considered prior to the implantation of the second bioprosthesis in order to guarantee the safety of the patient. A valve-in-valve should generally be avoided in patients with any of the following conditions:

- Degenerated bioprosthesis with a significant paravalvular leak (between vessel wall and prosthesis), not securely fixed in the vena cava, or is not structurally intact.
- Degenerated bioprosthesis with a partially detached leaflet that may break loose or obstruct the device function.

4.5. PRE-IMPLANTATION REQUIREMENTS

Transfemoral caval valve implantation shall be performed in a catheter laboratory or hybrid room with hemodynamic monitoring, high-quality fluoroscopic and echocardiographic imaging capabilities.

The procedure for the transfemoral valve implantation consists of a standard catheterization of the right or left femoral vein followed by the introduction of the delivery system in the vasculature for retrograde access to the vena cava superior or vena cava inferior and deploying the bioprosthesis in the vena cava.

It is necessary to perform venous angiography before insertion of the delivery system, to facilitate the placement of the bioprosthesis in its proper position.

4.6. MEASUREMENT PROTOCOL

The selection of the correct bioprosthesis size requires the analysis of the measurements of the location of valve implants in the vena cava, for that it is necessary to use the following measurements descriptions. The images must be created using a Gated-ECG Cardiac Multidetector Computed Tomography (MDCT) or Cardiac Multi-slice Computed Tomography. The images must be transferred to a validated 3D Multi-Planar Reconstruction software (3D MPR), like Osiris, Horus, CVI42, 3Dmension, Terarecon and other similar programs.

CAUTION: Correct sizing of the bioprosthesis is essential to mitigate the risk of paravalvular leakage and migration of the valve.


4.6.1. MEASUREMENT OF THE SUPERIOR VENA CAVA (SVC)

A. MEASUREMENTS

For implantation of SVC bioprosthesis 7 measurements (all measurements in mm) are required as directed in the following checklist and Figure 3.

| No | Measurement | Measured value in mm |
|----|--|----------------------|
| 1 | Diameter of confluence | |
| 2 | Diameter of SVC at level of top of Pulmonary Artery | |
| 3 | Diameter of SVC at level of middle of Pulmonary Artery | |
| 4 | Diameter of SVC at level of bottom of Pulmonary Artery | |
| 5 | Diameter of SVC-Right Atrium junction | |
| 6 | Length between measured point 1 and point 3 | |
| 7 | Length between measured point 1 and point 5 | |

Figure 3: SVC measurements

B. CALCULATE THE DIAMETER VARIATION

- Between 2 and 3;
- Between 3 and 4.

| No | Measurement | Value in mm |
|----|---|-------------|
| 1 | Difference in Diameter at measuring point 2 and point 3 | |
| 2 | Difference in Diameter at measuring point 3 and point 4 | |

C. HOW TO EXECUTE THE MEASUREMENT

Using the MPR plane, find the best projection at interest point of measurement and adjust the planes to be very orthogonal to SVC.

1. Put the center of line's intersection at the interest point. Repeat for the 3 planes presented in 3D MPR screen (Figure 4).
2. Adjust the coronal and sagittal view, rotating the axis of image to obtain an orthogonal (perpendicular) view in these both visualizations. At the end, use the adjusted axial view to measure the perimeter and area of SVC in different interest points (use 10 to 15 points to create the polygon for area and perimeter measurement).


Figure 4: Adjust planes

3. Use the software to measure the distance between different planes of measurement. If the SVC presents angulation or tortuosity, use the open polygon tool to trace the length at the center of vein. The lengths between confluence and point of SVC at middle of pulmonary artery and between confluence and SVC-RA junction are fundamental to implantation planning (Figure 5).


Figure 5: Trace length at center of vein

D. DIAMETER VARIATION

Use the measurements of SVC at three PA levels (top, middle and bottom) and distance between them (Figure 6).


Figure 6: Diameter variation

E. CALCULATE THE RELATION

- Diameter of SVC@TOP PA – Diameter of SVC@Middle PA)/Length between TOP PA and Middle PA.
- Diameter of SVC@Middle PA – Diameter of SVC@Bottom PA)/Length between Middle PA and Bottom PA.

The results from this calculation will be the rate of dilation of SVC per mm.

F. HOW TO CHOOSE THE SVC VALVE

To select the valve for SVC to be implanted it is obligatory to observe the valve measurements (Figure 7) and the Table 2 with the limits of the measurements.


Figure 7: TricValve® for SVC - measurements

| Measurement | SVC 25 | SVC 29 |
|----------------------|---------------------|---------|
| Confluence | Larger than 14mm | |
| SVC TOP PA | 19-31mm | 22-34mm |
| SVC Middle PA | 22-31mm | 27-34mm |
| Diameter variation 1 | Less than 0.50mm/mm | |
| Diameter variation 2 | | |
| Length Middle PA | Larger than 35mm | |
| Length to SVC-RA | Larger than 50mm | |

Table 2: TricValve® for SVC - measurements

4.6.2. MEASUREMENT OF THE INFERIOR VENA CAVA (IVC)

A. MEASUREMENTS

For implantation of the IVC valve 5 measurements (all measurements in mm) are required as directed in the following checklist and Figure 8.

1. IVC-RA transition diameter;
2. IVC at top of Hepatic Veins;
3. IVC just below Hepatic Veins;
4. IVC at 5cm below IVC-RA transition;
5. IVC-RA transition (1) to IVC top of HV (2) length.


Figure 8: IVC measurements

| No | Measurement | Measured Value in mm |
|----|---|----------------------|
| 1 | IVC-RA transition diameter | |
| 2 | IVC at top of Hepatic Veins | |
| 3 | IVC just below Hepatic Veins | |
| 4 | IVC at 5cm below IVC-RA transition | |
| 5 | IVC-RA transition (1) to IVC top of HV (2) length | |

B. HOW TO EXECUTE THE MEASUREMENT

Using the MPR plane, find the best projection at interest point of measurement and adjust the planes to be very orthogonal to IVC.

1. Put the center of line's intersection at the interest point. Repeat for the 3 planes presented in 3D MPR screen (Figure 9).
2. Adjust the coronal and sagittal view, rotating the axis of image to obtain an orthogonal (perpendicular) view in these both visualizations. At the end, use the adjusted axial view to measure the perimeter and area of IVC in different interest points (use 10 to 15 points to create the polygon for area and perimeter measurement).


Figure 9: Adjust planes

3. Use the software to measure the distance between different planes of measurement. If the IVC presents angulation or tortuosity, use the open polygon tool to trace the length at the center of vein. The lengths between junction of IVC to RA and top of hepatic veins is fundamental to implantation planning (Figure 10).


Figure 10: Trace length at center of vein

C. HOW TO CHOOSE THE TRICVALVE® FOR IVC

To select the valve for IVC to be implanted it is obligatory to observe the valve measurements (Figure 11) and the Table 3 with the limits of the measurements.


Figure 11: TricValve® for IVC - measurements

| Measurement | IVC 31 | IVC 35 |
|----------------------------------|------------------|------------|
| IVC-RA junction | 24 to 31mm | 28 to 35mm |
| IVC top HV | 24 to 31mm | 28 to 35mm |
| Length IVC/RA junction – H Veins | Larger than 10mm | |
| IVC just below HV | 21 to 35mm | 27 to 43mm |
| IVC 5cm below RA junction | 21 to 35mm | 27 to 43mm |

Table 3: TricValve® for IVC - Measurements

5. REQUIRED PRODUCTS

- TricValve® Transcatheter Bicaval Valves System for SVC
- TricValve® Transcatheter Bicaval Valves System for IVC

The TricValve® Transcatheter Bicaval Valves System is packaged in a blister, sealed in Tyvek pouch configuration and sterilized with ethylene oxide (EO). The bioprosthesis is supplied already loaded inside the delivery system. The delivery system is in sterile condition if the pouches are not damaged or opened. The outer surfaces of the outer bag are NOT STERILIZED and should not be placed in the sterile field.

6. EXPERIENCE REQUIRED BY PHYSICIAN

Physicians using the Bioprosthesis and delivery system must have experience in

- Caval valve implantation (CAVI)
- Transfemoral access and catheterization
- Self-expandable vena cava valve bioprosthesis implantation using transcatheter/transfemoral procedures.

The implantation of TricValve® Transcatheter Bicaval Valves should only be performed by doctors who have received a P&F Products & Features TricValve® Training.

7. ADDITIONAL MATERIALS

The materials necessary for the implantation of the Bioprosthesis are:

- standard cardiac catheterization lab equipment;
- sterile syringe, 20 ml;
- Fluoroscopy (fixed, mobile or semi-mobile fluoroscopy systems appropriate for use in percutaneous coronary interventions)
- Transesophageal (TEE) or Transthoracic Echocardiography capabilities;
- 2x Sheath 6F;
- 1x Standard Guidewire 0.035”;
- 1x Amplatz Super stiff guidewire;
- 1x Dilator 14 F;
- 1x 5F Pigtail;
- 1x Standard right-heart balloon-Catheter (as used for right heart Cath);
- 2x Pro-glide (optional).

This material is not part of the TricValve® Transcatheter Bicaval Valves System and must be supplied by the hospital prior to the surgical procedure.

8. INSTRUCTIONS FOR USE

CAUTION: Once the TricValve® Transcatheter Bicaval Valves System is removed from the packaging, ensure all subsequent procedures are performed in a sterile field.

CAUTION: Correct sizing of the bioprosthesis is essential to mitigate the risk of paravalvular leakage and migration of the Bioprosthesis.

8.1. INSPECTION PRIOR TO USE

1. Before removing the TricValve® Transcatheter Bicaval Valves System from its primary packaging, carefully inspect the packaging for any evidence of damage that could compromise the sterility or integrity of the device (e.g., broken or missing seals, torn or punctured pouch).

CAUTION: Do not use the product after the expiry date has elapsed or if the integrity of the sterilized packaging has been compromised (e.g. damaged packaging).

2. Visually check that the product is free of defects. Do not use if any defects are noted.


8.2. PREPARATION OF THE TRICVALVE® TRANSCATHETER BICAVAL VALVES SYSTEM

Wipe the length of the catheter with a moist (saline) gauze.

1. The catheter with the loaded bioprosthesis is slightly open. Use the micro knob on the handle to close the catheter and then open until the tube extremity reaches the edge of the bioprosthesis.

2. Attach a stopcock to the first flush port. Attach a 20-mL syringe filled with saline to the stopcock on the first flush port and flush. Repeat step until removing air bubbles from the tube.

Figure 12:
Attach stopcock
and flush

NOTE: The bioprosthesis leaflets are dehydrated when loaded in the delivery system and hence requires hydration using normal saline to achieve its required state for implantation. The hydration procedure is done by flushing the device with normal saline as per above mentioned procedure. It's important to remove all air bubble and flush with minimum of 60mL of saline solution and keep the bioprosthesis wet during minimum 2 minutes.

3. Verify no catheter leakage is observed during any of the flushing steps. If leakage is found, use a new system.

4. Attach a 10-mL syringe filled with saline to the guidewire port on the handle on the catheter and flush.

5. Close the capsule using the micro knob on the backside of the handle until the extension of the bioprosthesis is covered.

6. If necessary, flush the first port to remove any residual air bubble.

7. Connect the flushing line to the pressure line.

8. Conduct a final visual inspection of the loaded bioprosthesis to make sure the frame is free of creases or in folds beyond the second node from the inflow end. Ensure check is performed circumferentially around the entire bioprosthesis.

CAUTION: If a crease or fold greater than 2 nodes long is noticed, do not use the bioprosthesis or delivery system.

Figure 13:
Close capsule



Figure 14:
Connect pressure line

8.3. VASCULAR ACCESS

NOTE: Vascular access should be achieved per hospital protocol (either percutaneously or via surgical cut down).

NOTE: The primary access vein will be used to introduce the valve and the delivery system; the secondary access vein will be used to introduce the reference pigtail.

1. Insertion of temporary pacemaker catheter if necessary.
2. Insert a 6-Fr introducer sheath into the secondary access vein.
3. A sheathless approach is also possible according to operator's distinction.
4. Administer anticoagulant according to hospital protocol. If heparin is administered as an anticoagulant, check the activated clotting time (ACT) after initial bolus of heparin and recheck every 30 minutes thereafter. Maintain ACT \geq 250 seconds.

NOTE: Anticoagulant may be administered at any time prior to this point, but avoid delaying beyond this point.

8.4. IMPLANTATION PROCEDURE

1. Place a 5F arterial catheter in the left femoral artery used for continuous blood pressure monitoring.
2. Position a pulmonary catheter through the left femoral vein in the right pulmonary artery (rPA) to mark the crossing of the rPA with the SVC.
3. Introduce a 6F pigtail catheter through a 6F sheath in the right femoral vein.
4. Obtain an angiogram of the SVC prior to bioprosthesis deployment.
5. Exchange the pigtail for a straight 0.035 in. stiff guidewire with a soft tip.
6. Make a small skin incision.
7. Advance the TricValve® Delivery System over the guide wire through the femoral vein and the Inferior Vena Cava (IVC) into the RA and Superior Vena Cava (SVC).
8. Place the upper part of the SVC bioprosthesis in the confluence, with the belly of the SVC bioprosthesis positioned above the rPA crossing.
9. Confirm the catheter position under fluoroscopic and echocardiographic visualization and unload the uppermost 20 mm of the valve partially.
10. After the position is confirmed under fluoroscopy, fully deploy the SVC bioprosthesis by using the micro knob.
11. Retrieve the TricValve® Delivery System with the guidewire kept in place.
12. Pressure measurements may be taken to ensure correct prosthetic valve function.
13. Withdraw the catheter in the rPA to avoid interference with the IVC valve.
14. Insert the TricValve® Delivery System with the already loaded IVC bioprosthesis at the puncture site.
15. Position the IVC bioprosthesis at the height of the diaphragm with the skirt visible just above the hepatic vein inflow.
16. Align the constrained segment of the stent frame with the cavo-atrial junction by careful pullback of the catheter.
17. With a safety margin of 5 mm, take care to avoid a low or high valve position causing either hepatic vein obstruction or paravalvular regurgitation.
18. After a satisfactory position is confirmed, release the IVC bioprosthesis slowly from the catheter, observing the distal deployment. This avoids a jerky release of the bioprosthesis.
19. Under fluoroscopic guidance, confirm that the catheter tip is coaxial with the inflow portion of the bioprosthesis.
20. Withdraw the catheter to the femoral vein while maintaining guidewire position.
21. Close the capsule and remove the catheter through the femoral access.

NOTE: If the capsule does not close properly, gently rotate the catheter clockwise ($<180^\circ$) and then counter clockwise ($<180^\circ$) until the capsule closes.

CAUTION: Ensure the capsule is closed before catheter removal. If increased resistance is encountered when removing the catheter through the introducer sheath, do not force passage. Increased resistance may indicate a problem and forced passage may result in damage to the device and/or harm to the patient. If the cause of resistance cannot be determined or corrected, remove the catheter and introducer sheath as a single unit over the guidewire, and inspect the catheter and confirm that it is complete.

8.5. POST DEPLOYMENT

1. Remove the reference pigtail catheter over a standard guidewire.
2. Remove the 6-Fr introducer and close the access site per hospital protocol.
3. Administer anticoagulation and/or antiplatelet therapy as required according to hospital protocol.
4. Perform Transesophageal echocardiography (TEE Echo) immediately after the procedure, if needed.

8. FINAL RESULT


Figure 15 - Implanted bioprosthesis

9. MEDICAL TEAM TRAINING

The P+F Products + Features GmbH provides detailed TricValve® Training to operators. For on-site training and assistance materials, contact a P+F Products + Features GmbH representative (see www.productsandfeatures.com).

CAUTION: The implantation of TricValve® Transcatheter Bicaval Valves System should only be performed by doctors who have received a P+F Products + Features GmbH TricValve® Training.

10. PACKAGING

The TricValve® Transcatheter Bicaval Valves System is supplied sterile and non-pyrogenic in individual packaging. The Bioprosthesis is already loaded inside the delivery system. The TricValve® Transcatheter Bicaval Valves System is packaged in a double pouch sealed Tyvek bag and sterilized with ethylene oxide.

WARNING: The outer surface of the sealed pouch is not sterile and therefore shall not come into contact with sterile field and instruments.

11. STORAGE

Store the TricValve® Transcatheter Bicaval Valves System in a dry area at room temperature (8°C to 25°C [46,4°F to 77°F]). Avoid exposing to extreme temperature fluctuations and direct contact with sunlight. Avoid freezing as it can damage the bioprosthesis. Storage at elevated temperatures may potentially damage the polymeric components and adhesives thus compromising the product performance.

12. MRI SAFETY

Non-clinical testing has demonstrated that the valve is MRI conditional. A patient with this device can be safely scanned in an

MR system meeting the following conditions:

- Static magnetic field of 1.5 tesla (T)
- Spatial gradient field of ≤ 2500 gauss/cm (25 T/m).
- Normal operating mode only with a maximum whole body specific absorption rate (SAR) of 2.0 W/kg for 15 minutes as read from equipment monitor.

In non-clinical testing, the image artifact caused by the device extends as far as 14.5 mm from the implant for spine echo images and 30 mm for gradient echo images when scanned in a 3.0 T MRI system. The artifact obscures the device lumen in gradient echo images. The implant has not been evaluated in MR systems other than 1.5 or 3.0 T.

13. PATIENT REGISTRATION

A patient registration form is included in each TricValve® Transcatheter Bicaval Valves System package. After implantation, all requested information must be completed. The serial number is located on the package. Return the original form to P+F Products + Features GmbH address indicated on the form and provide the temporary identification card to the patient prior to discharge.

P+F Products + Features GmbH will provide an Implanted Device Identification Card to the patient. The card contains the name and telephone number of the patient's physician as well as information that medical personnel would require in the event of an emergency.

14. DISPOSAL OF USED CLINICAL DELIVERY SYSTEMS

The used delivery system may be disposed in the same manner that the hospital waste and biohazardous materials are handled. There are no special or unusual risks related to the disposal of the delivery system which are required to be followed.

15. RETURN OF EXPLANTED VALVES

P+F Products + Features GmbH is interested in obtaining explanted TricValve® Transcatheter Bicaval Valves. The explanted bioprosthesis should be placed into a suitable histological fixative such as 10% formalin or 2% glutaraldehyde and returned to P+F Products + Features GmbH. Refrigeration is not necessary under these circumstances. Specific pathological studies of the explanted valve will be conducted under the direction of a specialized pathologist. A written summary of the findings will be returned to the physician. For a product return kit, contact a P+F Products + Features GmbH distribution center or a representative. For further instructions on the return of an explanted device, please contact P+F Products + Features GmbH.

16. DISCLAIMER OF WARRANTY

ALTHOUGH THE P+F PRODUCTS + FEATURES GMBH' TRICVALVE® BICAVAL TRANSCATHETER VALVE SYSTEM, HEREAFTER REFERRED TO AS "PRODUCT", HAS BEEN CAREFULLY DESIGNED, MANUFACTURED, AND TESTED PRIOR TO SALE, THE PRODUCT MAY FAIL TO PERFORM IT'S INTENDED FUNCTION SATISFACTORILY FOR A VARIETY OF REASONS. THE WARNINGS CONTAINED IN THE PRODUCT LABELING PROVIDE MORE DETAILED INFORMATION AND ARE CONSIDERED AN INTEGRAL PART OF THIS DISCLAIMER OF WARRANTY. P+F PRODUCTS + FEATURES GMBH, THEREFORE, DISCLAIMS ALL WARRANTIES, BOTH EXPRESS AND IMPLIED, WITH RESPECT TO THE PRODUCT. P+F PRODUCTS + FEATURES GMBH SHALL NOT BE LIABLE FOR ANY INCIDENTAL OR CONSEQUENTIAL DAMAGES CAUSED BY ANY USE, DEFECT OR FAILURE OF THE PRODUCT, WHETHER THE CLAIM IS BASED ON WARRANTY, CONTRACT, TORT OR OTHERWISE.

The exclusions and limitations set out above are not intended to, and should not be construed so as to contravene mandatory provisions of applicable law. If any part or term of this DISCLAIMER OF WARRANTY is held to be illegal, unenforceable or in conflict with applicable law by a court of competent jurisdiction, the validity of the remaining portions of this DISCLAIMER OF WARRANTY shall not be affected, and all rights and obligations shall be construed and enforced as if this DISCLAIMER OF WARRANTY did not contain the particular part or term held to be invalid.

17. LABELING

| | |
|---|---------------------|
| REF | Reference number |
| LOT | Lot number |
| SN | Serial number |
|  | Date of manufacture |

| | |
|---|--|
|  | Use by Date |
|  | Manufacturer |
|  | Consult Instructions for Use |
|  | Caution |
|  | Size |
|  | Minimum Sheath Size |
|  | Recommended Guidewire Size |
|  | Do not use if package is damaged |
|  | Keep away from Sunlight |
|  | Keep dry |
|  | Temperature limit |
|  | Sterilized with ethylene oxide |
|  | Do not resterilize |
|  | Single Use, Do not reuse |
|  | Non pyrogenic |
|  | MR conditional |
|  | Maximum Guidewire Diameter |
|  | Quantity |
|  | Conformité Européenne (European Conformity). This symbol means that the device fully complies with applicable European Union Acts. |

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© 2019 P&F Products & Features Germany | All rights reserved.