
INSTRUCTIONS FOR USE

TricValve® Transcatheter Bicaval Valve System (SVC 25, SVC 29, IVC 31, IVC 35)

ENGLISH	EN	Instructions for Use
GERMAN	DE	Gebrauchsanweisung
FRENCH	FR	Mode d'emploi
SPANISH	ES	Instrucciones de uso
ITALIAN	IT	Istruzioni per l'impiego
DUTCH	NL	Instructies voor gebruik
POLISH	PL	Sposób stosowania
PORTUGUESE	PT	Instruções de Utilização
SWEDISH	SV	Instruktioner för användning
FINNISH	FI	Käyttöohjeet
ROMANIAN	RO	Instrucțiuni de utilizare
CZECH	CS	Návod k použití
HUNGARIAN	HU	Használati útmutató
BULGARIAN	BG	Инструкции за употреба
GREEK	EL	Οδηγίες Χρήσης
SLOVAK	SK	Návod na použitie
LITHUANIAN	LT	Naudojimo instrukcijos
NORWEGIAN	NO	Bruksanvisning
TURKISH	TR	Kullanım Kılavuzu
RUSSIAN	RU	Инструкции по применению
SIMPLIFIED CHINESE	ZH	使用说明
AZERBAIJANI	AZ	İstifadə üçün təlimatlar

CAUTION: Carefully read these instructions for use, which address the warnings, precautions, and residual risks for this medical device. Implantation of TricValve® Transcatheter Bicaval Valve System should only be performed by physicians who have received or undergone P+F Products + Features GmbH training on device usage.

- Do not attempt to use the TricValve® Transcatheter Bicaval Valve (herein referred to as the TricValve®) or TricValve® Transcatheter Bicaval Valve delivery system (herein referred to as TricValve® delivery system) before completely reading and understanding the information contained in this booklet.
- Failure to follow these instructions, warnings, and precautions may lead to device damage or patient injury. Use of the TricValve® system should be restricted to those physicians trained to perform invasive endovascular procedures and to those physicians trained in the proper use of the system.
- Consult with authorized P+F Products + Features Technical personnel for appropriate selection of TricValve® valve size.
- The TricValve® Transcatheter Bicaval Valve System (IVC and SVC) are supplied STERILE using EO gas sterilization process. The outer Tyvek® pouch used to package the device is NON-STERILE.
- All TricValve® devices are supplied for single use only. After use, dispose in accordance with hospital administrative and/or government policy. Do not re-sterilize.

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

1 INTRODUCTION

The TricValve® Transcatheter Bicaval Valve System (herein referred to as the TricValve® system) is designed to be used in patients with symptomatic, severe tricuspid valve regurgitation without the need for conventional open-heart surgery. The TricValve® system consists of the following four (4) models – 2 for SVC implantation (SVC25 and SVC29) and 2 for IVC implantation (IVC31 and IVC 35), as outlined in below Table 1:

Table 1 Model / Sizes of TricValve® Transcatheter Bicaval Valve - REF Numbers

Device name / Model Name	REF Number
TricValve® Transcatheter Bicaval Valve System - SVC 25	PF01020025
TricValve® Transcatheter Bicaval Valve System - SVC 29	PF01020029
TricValve® Transcatheter Bicaval Valve System - IVC 31	PF01020031
TricValve® Transcatheter Bicaval Valve System - IVC 35	PF01020035

Caution: Safety and effectiveness for implantation of a single valve and for more than two valves have not been established.

2 MEDICAL DEVICE DESCRIPTION

The TricValve® Transcatheter Bicaval Valve system consists of 2 bioprostheses (TricValve® IVC and TricValve® SVC Valves) pre-loaded on two identical TricValve® delivery systems. The Medical Device is a bicaval transcatheter tricuspid valve implantation system, which includes the TricValve® Transcatheter Bicaval Valve for Superior vena cava (SVC) and for Inferior vena cava (IVC). The TricValve® Transcatheter Bicaval Valves (IVC and SVC) are individually packaged into two separate boxes, sterile and ready to use. The TricValve® system percutaneously delivers the implant to the defined location via femoral vein access using a transvenous (tricuspid) approach.

2.1. TRICVALVE® (BIOPROSTHESES)

The TricValve® bioprosthesis is developed for treatment of caval reflux present in cases of severe tricuspid regurgitation with a high risk for open-heart surgery. The TricValve® bioprosthesis is designed to treat severe tricuspid regurgitation without removal / disturbing the defective tricuspid valve. The TricValve® bioprosthesis is deployed in Inferior vena cava and Superior vena cava to treat regurgitation of tricuspid valve. The TricValve® bioprosthesis consists of a tubular metallic structure of nitinol which is self-expandable and radiopaque with three valve leaflets and skirt of bovine pericardium sutured and complemented by an additional skirt of polyester (IVC Valve only) to avoid paravalvular leaks. The bioprostheses leaflets are processed with anti-calcification as well as chemical dehydration treatment. The bioprosthesis is sterilized by Ethylene Oxide (EO) gas sterilization process and pre-loaded on the delivery system. The Shape of TricValve® bioprosthesis is designed to allow safe valve deployments in the vena cavae.

The TricValve® bioprostheses are available in different diameters for each model (SVC and IVC) specifically designed to adapt to the anatomic features of the superior and Inferior vena cava (Figure 1 and Table 2). The valve is available in different diameters as shown in Table 2. Valve size recommendations are based on the measurements of IVC and SVC, as measured by computed tomography (CT scan). Patient anatomical factors and imaging modalities

should be considered during the TricValve® size selection. Refer to the TricValve® physician training material for additional information on the selection of size of bioprostheses, including additional considerations for implant selection.

Table 2 Sizes of TricValve® Transcatheter Bicaval Valves (mm)

TRICVALVE® MODEL	VALVE SIZE (mm)	PROXIMAL DIAMETER (mm)	DISTAL DIAMETER (mm)	NOMINAL HEIGHT (mm)	VENA CAVA DIAMETER RANGE (mm)
SVC 25	25	25	20	67	22 – 31*
SVC 29	29	29	20	70	27 – 34*
IVC 31	31	34	38	66	24 – 31**
IVC 35	35	38	47	66	28 – 35**

* - Measured at pulmonary artery level, ** - Measured at hepatic vein level

Note: TricValve® Transcatheter Bicaval Valve System comprises TricValve IVC and SVC bioprostheses each preloaded on a TricValve® delivery system.

2.2. TRICVALVE® DELIVERY SYSTEM (CATHETER)

The TricValve® Transcatheter Bicaval Valve delivery system (Figure 2 and 3) is an enclosure retraction type system (28 Fr) and used to deploy the TricValve® bioprostheses. TricValve® bioprosthesis has anchors for mounting into delivery system for safe and accurate deployment at desired implant location. The proximal end of the system has an atraumatic radiopaque tip and a protective sheath. The protective sheath covers and maintains the valve in a crimped position. A stabilizer tube is fixed at the handle and extends outside the delivery system shaft. It provides a barrier between the inner delivery system shaft and vessel walls, thus enabling the delivery system to retract freely. The TricValve® delivery system is compatible with a 0.035" (0.889mm) guide wire.

The TricValve® delivery system handle includes a macro slider for opening and closing the valve capsule and a micro knob facilitate the accurate release of the TricValve® bioprosthesis. The micro knob rotates clockwise to open the valve capsule and in anticlockwise direction to close the valve capsule. The delivery system has a tube flush port which is used to hydrate the bioprostheses leaflets and remove air before usage.

There is a single TricValve® delivery system model which is used for delivery of both valve models.

3 FACILITIES AND TRAINING REQUIREMENTS

Facilities that intend to perform an implantation procedure utilizing the TricValve® Transcatheter Bicaval Valve system must have access to cine, fluoroscopy and transesophageal echocardiography (TEE) imaging capabilities throughout the procedure. In addition, the implanting physicians must have prompt access to facilities with the necessary equipment, instruments, supplies, and personnel to perform emergency tricuspid valve surgery, if required. A comprehensive training program is established by P+F Products + Features GmbH and must be completed by the implanting physicians before using the TricValve® system. The implanting physicians should have advanced technical knowledge and experience in related catheter-based procedures.

Physicians using the bioprostheses and delivery system must have experience in:

- Structural Heart Treatment
- Transfemoral access and catheterization
- Self-expandable vena cava valve using transcatheter/transfemoral procedures.

The implantation of TricValve® Transcatheter Bicaval Valves should only be performed by doctors who have received P+F Products + Features TricValve® Training.

4 INDICATIONS OF USE

4.1. INTENDED USE

The TricValve® Transcatheter Bicaval Valve System is used for the implantation of two dedicated self-expanding biological valves to treat symptoms of right heart failure, in patients with severe tricuspid regurgitation and significant caval reflux.

The TricValve® delivery system percutaneously delivers the TricValve® bioprosthesis via femoral vein access using a transvenous (tricuspid) approach.

4.2. INDICATIONS OF USE

The TricValve® Transcatheter Bicaval Valve System is indicated for patients with symptomatic severe tricuspid regurgitation (hemodynamically relevant) and caval reflux, in patients at high risk for open surgical therapy. The Medical Device is indicated for patients, as determined by a multidisciplinary team of cardiologists experienced in the evaluation and treatment of heart failure and tricuspid valve disease.

Since the device does not interfere with the natural tricuspid valve, it can also be implanted specifically in patients deemed surgically inoperable by the team of cardiologist. Patients with previous valve interventions (surgical or minimally invasive) or patients who are implanted with pacemaker are also eligible for TricValve® implantation.

4.3. TARGET PATIENT POPULATION

The risks and benefits of TricValve® Transcatheter Bicaval Valve System should be carefully considered for each patient by a team of cardiologists before use. The risks of long-term anticoagulant therapy should also be taken into account. For details about the intended use, indications, and contraindications, refer to the specific sections. TricValve® Transcatheter Bicaval Valve System is indicated for patients with severe tricuspid regurgitation (assessed with echo) and significant caval reflux.

TricValve® is intended for patients:

- With significant backflow in venous system (Right heart catheter V-wave $\geq$ 15 mmHg),
- Without precapillary pulmonary hypertension (Systolic pulmonary artery pressure sPAP <65 mmHg.)

- With remaining right ventricular function (Tricuspid annular plane systolic excursion TAPSE ≥ 12 mm),
- With Left Ventricular Ejection Fraction (LVEF) $\geq 30\%$,
- Not deemed suitable for cardiac surgery, including minimally invasive cardiac surgery,
- With an anatomy suitable for the available TricValve sizes, as assessed by CT or MRI

TricValve® patients' eligibility is not limited by the tricuspid annulus size, morphology or right atrium/right ventricle size and morphology. The presence of pacemaker leads is not a contraindication for the treatment with TricValve® Transcatheter Bicaval Valve System.

5 CONTRAINDICATIONS

The TricValve® Transcatheter Bicaval Valve system is contraindicated in patients who have any of the following conditions:

- Active gastritis or peptic ulcer or any other condition with an increased bleeding risk.
- Active endocarditis within 3 months or infection requiring antibiotic therapy (oral or intravenous) within 2 weeks of the scheduled implant.
- Evolutionary or recent stroke / Cerebrovascular accident (CVA).
- Creatinine clearance < 20 ml / min.
- Bleeding diathesis or coagulopathy or patients refusing blood transfusion.
- Untreatable hypersensitivity, allergy or contraindication to device's components, (including nitinol, nickel, titanium, bovine pericardium, polyester, PTFE, silicon, and/or polymeric materials).
- Patients in whom transesophageal echo (TEE) and transthoracic echo (TTE) is contraindicated during the procedure or screening TEE and TTE is unsuccessful.
- Patients' anatomy not suitable for available TricValve® sizes.
- Recent myocardial infarction (< 30 days)
- Sepsis, including active endocarditis.
- Thrombosis of the lower venous system or vena cava filter.
- Tricuspid valve anatomy that precludes proper device deployment and functionality based on CT and echocardiographic evaluation.
- Untreatable hypersensitivity or contraindications to any of the following: Vitamin K antagonists, all antiplatelets, heparin and all anticoagulants, or sensitivity to contrast media that cannot be adequately premedicated.
- Vascular conditions (e.g. stenosis, tortuosity) that make insertion and endovascular access impossible to the superior or Inferior vena cava.

Note: Patient must be able to tolerate at least one antiplatelet medication AND one anticoagulant medication.

6 WARNINGS

Note: Before use, carefully read all warnings, and instructions for use for all components of the system. Failure to read and follow all instructions or failure to observe any stated warnings/instructions could cause serious injury or death to the patient.

- Do not use TricValve® Transcatheter Bicaval Valve System outside of the labelled indication.

- The TricValve® devices are designed, intended, and distributed for single use only. Do not re-sterilize or re-use the devices. There is no data to support the sterility, non- pyrogenic, and functionality of the devices after reprocessing or re-sterilization. Such action could lead to illness or an adverse event as the device may not function as originally intended.
- Ensure proper sterile techniques are utilized during the preparation, transfer, and use of the devices.
- Do not use the system if the tamper evident seal is broken or the valve is damaged, or the expiration date has elapsed.
- Do not use the system if the "Use by date" or "Expiry Date" has passed.
- Do not use the catheter if the packaging or any components are not sterile, which may not be apparent under visual inspection, may result. This will help prevent leaflet damage that may impact valve functionality.
- Caution should be exercised in removing the system from the package, to prevent any damage.
- TricValve® cannot be exposed to solutions, antibiotics, or chemicals other than sterile physiological saline solutions. Irreparable damage to the device, which may not be apparent under visual inspection, may result. This will help prevent leaflet damage that may impact valve functionality.
- Accelerated deterioration of the TricValve® bioprostheses may occur in patients with an altered calcium metabolism.
- Verify that relevant patient anatomical parameters are suitable and within the specifications for performing the procedure and implantation of TricValve® bioprostheses.
- Ensure correct valve sizing by evaluation of CT images with P+F technical team. Incorrect sizing of the bioprostheses may lead to paravalvular leak (PVL), device migration, embolization and/or cardio-vascular rupture.
- Do not use it if any of the devices were dropped, damaged, or mishandled during any part of the procedure. This will require replacement of the TricValve® System.
- Do not handle or manipulate the TricValve® bioprostheses with sharp or pointed objects.
- Do not open the contents of the package until you are sure about the deployment and the appropriate size of the bioprosthesis. This can compromise the sterile conditions of medical devices.
- Standard flushing and de-airing techniques should be used during preparation and throughout procedures to prevent air embolism.
- After the delivery system was inserted into patient and release of valve is unsuccessful / out of capsule; do not attempt to reload the same. Dispose the entire system and do not re-use any component.
- Patients with pre-existing cardiac leads should be carefully assessed prior to implantation to avoid potential adverse interaction between devices. Care should be taken when implanting cardiac leads after TricValve® implantation to avoid potential adverse interaction between the devices.
- Patients with previously implanted devices (e.g., IVC filter) should be carefully assessed prior to insertion of the delivery system to avoid potential damage to vasculature or a previously implanted device.
- Caution should be exercised during implanting a valve in patients with clinically significant coronary artery disease. Caution should be exercised in patients with relevant left side valve disease.
- Special caution should be exerted In patients with patent foramen ovale and atrial septal defect with significant shunt, closure of the defect should be considered before TricValve implantation.

- As with any implanted medical device, there is a potential for an adverse immunological response.
- Serious adverse events, sometimes leading to surgical intervention and/or death, may be associated with the use of this system ("Potential Adverse Events"). A full explanation of the benefits and risks should be given to each prospective patient before use.
- The presence of the following conditions may be associated with higher procedural morbidity:
 - Unstable angina
 - Active infection
- Careful and continuous medical follow-up is advised so that implant-related complications can be diagnosed and properly managed.
- Patients implanted with the TricValve® should maintain anticoagulant/antiplatelet therapy as determined by their physicians in accordance with current guidelines, to minimize the risk of valve thrombosis or thromboembolic events.
- The TricValve® Transcatheter Bicaval Valve System has not been evaluated in pregnant or pediatric patients (under 18 years of age).
- If any valve leaflets are mishandled or damaged during the procedure, replace the device with a fresh (new) one.
- Any mechanical failures of the delivery system may result in complications and severe injury to the patient.
- Do not use gloves with powder, to avoid any contamination of the TricValve® System.

7 PRECAUTIONS

Do not use the TricValve® Transcatheter Bicaval Valve System:

- If the sterile package is damaged.
- If the bioprosthesis is damaged.
- If the delivery system is damaged.
- If the delivery system is unable to flush or shows any leakage.
- If the guidewire does not slide easily inside the delivery system.
- If the use-by-date has elapsed.

7.1. PRECAUTIONS PRIOR TO USE

- The patient's eligibility for treatment with TricValve® device depends on the anatomic conditions based on CT scan.
- It is advised that a multi-disciplinary heart team carefully assesses eligibility and takes into consideration patient risk, anatomical suitability, and the severity of disease. The team shall assess the selected valve size for treatment and be of the opinion TricValve® implantation is preferable to alternative to percutaneous device solutions, including minimally invasive open-heart surgery.
- The reference catheter should be placed in the superior trunk of the right pulmonary artery (RPA) for positioning reference of the TricValve® Valve for the Superior vena cava (SVC).

7.2. REPOSITIONING PRECAUTIONS

- After initiation of the release, withdrawal of the bioprosthesis (e.g. catheter use) is not recommended. A withdrawal of a partially released bioprosthesis

- can cause mechanical failure of the releasing device, damage to a vena cava, an iliac or femoral vein, atrium, as well as vascular complications, prosthetic valve dysfunction (including an abnormal position of the device), embolization, and / or lead to emergency surgery.
- During release, the bioprosthesis may be retracted, if the vena cava has not been touched. After contact with the vein, it will not be possible to retract the bioprosthesis.
- After release, repositioning of the bioprosthesis is not recommended. Repositioning a released bioprosthesis may damage the vena cava or the atrium, as well as lead to vascular complications, prosthetic valve dysfunction (including an abnormal position of the device), embolization, and / or lead to emergency surgery.
- Do not try to remove the delivery system if some piece of the bioprosthesis is still connected with the delivery system. Ensure that there was a full bioprosthesis release prior to initiating the delivery system withdrawal procedure.

7.3. PRECAUTIONS AFTER IMPLANTATION

- After the valve is fully deployed, retrieve the delivery system tip by unlocking the micro knob, then, using your hand, hold the macro slider and slowly pull back the back-handle of the delivery system until the tip is fully retracted, then lock the micro knob. Make sure the valve capsule is closed before removing the delivery system. Keeping the guidewire in place, safely remove the TricValve® delivery system ensuring the tip is retracted.
- If you feel an increase in resistance when removing the delivery system, do not force the passage. Increased resistance may be indicative of a problem and may damage the device and / or cause injury to the patient if the passage is forced. If the cause of the resistance cannot be determined or corrected, slowly rotate the delivery system, and inspect whether the delivery system is complete or if there is any bent or ruptured region. One must always make sure that there is no connection of the delivery system with the valve during the removal of it.
- Care should be taken during removal of the delivery system from the patient, maintaining radioscopic visualization.
- Care should be taken when crossing a TricValve® Transcatheter Bicaval Valve implanted with ancillary devices.

CAUTION: To avoid damage to the device and/or injury to the patient, the process of tip retraction must be followed in its entirety and under fluoroscopy imaging guidance.

8 POTENTIAL ADVERSE EVENTS

Complications associated with standard cardiac catheterization, the use of anesthesia, and use of the TricValve® Transcatheter Bicaval Valve System could lead to the following outcomes: conversion to open surgery, emergent or non-emergent reoperation, explant, permanent disability, or death. Potential adverse events associated may include, but are not limited to, the following:

- Abnormal laboratory values (including electrolyte imbalance)
- Allergic reaction to anesthesia, contrast media, anticoagulation medications, or device materials

- Anemia or decreased Hemoglobin (Hgb) that may require transfusion.
- Aneurysm or pseudoaneurysm
- Angina or chest pain
- Arterio-venous fistula
- Asymptomatic cerebral infarct
- Bleeding
- Cardiac arrest
- Cardiac arrhythmia atrial (i.e., atrial fibrillation, supraventricular tachycardia)
- Cardiac arrhythmia ventricular (i.e., ventricular tachycardia, ventricular fibrillation)
- Cardiac insufficiency
- Cardiac tamponade / Pericardial effusion
- Cardiogenic shock
- Cardiovascular injury (including rupture, perforation or dissection of vessels, ventricles, atrium, or cardiac structures that may require intervention)
- Coagulopathy, coagulation disorder, bleeding diathesis
- Complications related to vascular access (e.g. dissection, perforation, pain, hemorrhage, hematoma, pseudo aneurysm, irreversible nerve damage, compartment syndrome, arteriovenous fistula, stenosis)
- Conduction system injury, which may require implantation of a pacemaker (temporary or permanent)
- Damage to or interference with function of pacemaker or implantable cardioverter defibrillator (ICD)
- Death
- Deep vein thrombosis (DVT)
- Deterioration of native valve (e.g. leaflet tearing, retraction, thickening)
- Diaphragmatic dysfunction
- Dislodgement of previously deployed implant
- Dyspnea
- Edema
- Electrolyte imbalance
- Emboli/embolization including air, particulate, calcific material, or thrombus.
- Embolization of the auxiliary device
- Encephalopathy
- Endocarditis
- Erosion of tissues
- Esophageal irritation
- Esophageal perforation or stricture
- Exercise intolerance or weakness
- Failure to retrieve any TricValve® system components.
- Fever
- Gastrointestinal bleeding or infarct
- General pain
- Hemorrhage requiring transfusion or intervention.
- Heart failure or low cardiac output
- Heart murmur
- Hematoma
- Hemodynamic compromise
- Hemolysis
- Hypersensitivity to device material like Nitinol, Bovine Pericardium, etc.
- Hypotension or hypertension
- Implant (valve) deterioration (wear, tear, fracture, or other)

- Implant malposition or failure to deliver to intended site.
- Infection (including septicemia and endocarditis)
- Inflammation
- Injury of the vena cava
- Intestinal ischemia
- Mesenteric ischemia or bowel infarction
- Migraine, transient ischemic attack (TIA) or other neurological deficits
- Multiple organ insufficiency
- Multi-system organ failure
- Myocardial infarction
- Myocardial ischemia
- Nausea and/or vomiting
- Necessity for the conversion to new non-emergency surgery
- Necessity for the conversion to non-emergent reoperation
- Necessity for the conversion to open heart surgery
- Necessity for the conversion to percutaneous valve intervention
- Neurological symptoms, including dyskinesia, without diagnosis of TIA or stroke.
- Neuropathic pain referred to shoulder subscapular pain.
- Pain (Shoulder and other body parts)
- Paralysis
- Paravalvular Leak (PVL)
- Percutaneous Valve intervention
- Pericardial effusion
- Peripheral ischemia
- Permanent disability
- Phrenic and other nerve injury
- Pleural effusion
- Pneumonia
- Prosthetic valve dysfunction including, but not limited to, fracture; bending (deformed conformation) of the bioprosthesis structure; sub expansion of the valve structure; calcification; pannus; wear; tear; prolapse or retraction in leaflets; poor valve coaptation; failure or disruption of suture; leakage; incorrect sizing (prosthesis-patient maladaptation); incorrect positioning (too high or too low) / incorrect placement; regurgitation and stenosis.
- Pulmonary Edema
- Pulmonary embolism
- Reaction to anti-platelet or anti-coagulation agents
- Renal insufficiency or renal failure (including acute renal injury)
- Respiratory compromise, respiratory insufficiency, or respiratory failure
- Retroperitoneal bleed
- Septal damage or perforation
- Septicemia, Sepsis
- Severe or small bleeding may or may not require transfusion or intervention (including fatal hemorrhage or disabling bleeding)
- Stroke
- Structural deterioration (wear, fracture, calcification, leaflet tear, leaflet thickening, stenosis of implanted device, or new leaflet motion disorder)
- Syncope
- Temporary phrenic nerve neuralgia
- Thromboembolism
- Thrombosis / clot (including valve thrombosis)
- Transient ischemic attack (TIA)

- Urinary tract infection and/or bleeding
- Valve dislodgement / embolization
- Valve endocarditis
- Valve injury
- Valve malposition
- Valve migration
- Valve stenosis
- Valve Thrombosis
- Valvular regurgitation
- Vascular injury or trauma, including dissection or occlusion.
- Ventricular wall damage or perforation
- Vessel spasm
- Worsening of heart failure
- Worsening regurgitation / valvular insufficiency
- Wound dehiscence delayed or incomplete healing.

For a patient/user/third party in the European Economic area; if, during the use of this device or as a result of its use, a serious incident has occurred, please report it to the manufacturer (vigilance@productsandfeatures.com) and your national competent authority, details for which can be found at: https://ec.europa.eu/growth/sectors/medical-devices/contacts_en.

9 ADDITIONAL EQUIPMENT

The implantation procedure requires additional equipment that is not supplied with the TricValve® Transcatheter Bicaval Valve System. The additional equipment is provided below.

9.1. EQUIPMENT / ACCESSORIES REQUIRED

- Standard cardiac catheterization lab equipment
- High-Quality fluoroscopy imaging system
- Transesophageal echocardiography (TEE) or Transthoracic echocardiography (TTE) equipment
- Hemodynamic monitoring equipment
- Femoral introducer sheath
- Venous puncture kit - puncture needle, sheath, and guidewire
- 1x Set of step-up dilators from 14 to 27 Fr
- 1x Standard length 0.035-inch (0.89 mm) extra-stiff guidewire
- 1x Standard length 0.035-inch (0.89 mm) guidewire
- 1x 5 Fr Right coronary arterial catheter and guidewire
- 1x Pulmonary catheter
- 1x 6 / 5 Fr Pigtail angiographic catheter
- 2x Pro-glide (optional)
- 2x Sheath 6 Fr
- 3-way Stopcock
- 5 Fr Multipurpose catheter
- 1000 ml (minimum) sterile physiological saline solution
- Luer lock syringe (≥ 20 cc)
- Sterile syringe (≥ 20 cc)
- Diluted radiopaque contrast medium (15:85 medium to saline dilution)

These materials are ancillary to the TricValve® Transcatheter Bicaval Valve System and must be arranged / supplied by the hospital prior to the treatment procedure.

10 INSTRUCTIONS FOR USE

TricValve® Transcatheter Valve sizing for the procedure is performed by authorized P+F Products + Features personnel. The CT images of patients shall be used to generate Multiplanar Reconstruction (MPR) and calculations shall be done from these images. The Device size selection will be based on this report generated by the authorized P+F personnel for the Patient.

10.1. INSPECTION PRIOR TO USE

Step	Procedure
1	Before removing the TricValve® Transcatheter Bicaval Valve System from its primary packaging, carefully inspect the packaging for any evidence of damage that could compromise the sterility or integrity of device (e.g. broken or missing seals, torn or punctured pouch).
2	Visually check that the product is free of defects.

WARNING: Do not mishandle the delivery system or use the delivery system and accessory devices if the packaging and/or sterile barriers and any components have been opened or damaged, or the expiration date has elapsed, as sterility and/or function may be compromised.

WARNING: Do not use the TricValve® system if the expiration date has elapsed, as either sterility or valve function may be compromised.

WARNING: Do not use the TricValve® system if the tamper evident seal is broken, as sterility may be compromised.

CAUTION: Do not use the product if any defects are noted.

CAUTION: The TricValve® Transcatheter Bicaval Valve System is packaged in a double-pouch configuration and sterilized with Ethylene Oxide (EO) gas. The system is sterile if the pouches are undamaged and unopened. The outer surface of the outer pouch is NOT STERILE and must not be placed in the sterile field.

CAUTION: Once the TricValve® Transcatheter Bicaval Valve System is removed from the outer packaging, ensure all subsequent procedures are performed in a sterile field.


CAUTION: Correct sizing of the valves is essential to mitigate the risk of paravalvular leakage and migration of the valve.


CAUTION: The valves are dehydrated when loaded into the delivery system. Rehydration process is essential to ensure the hydration of the valve for its proper functionality.

10.2. DEVICE PREPARATION

All the device preparation steps are required to be completed for both IVC and SVC systems. It must be completed prior to implantation of TricValve® Transcatheter Bicaval Valve System. An implant card is provided in the main box for the patient.

Step	Procedure
1.	Carefully open the box package and remove the TricValve® package.
2.	Place the packed TricValve® system on table and inspect for any tears or damage to the outer non-sterile Tyvek® pouch.
3.	Carefully open the outer pouch and present the inner sterile Tyvek® package with TricValve® system to the scrub nurse / technician in the sterile field. The inner pouch should be handled using standard aseptic techniques to prevent contamination.
4.	The scrub nurse / technician in the sterile field shall open the inner package and carefully remove the blister tray containing the TricValve® Transcatheter Bicaval Valve System. Do not open the TricValve® inner pouch until implantation is certain. Caution: Avoid touching of the sterile contents in blister tray with non-sterile surfaces.
5.	Inspect TricValve® Transcatheter Bicaval Valve System for any signs of damage. The device should not be used for the treatment procedure and replaced with a fresh one if any damage is observed.
6.	Verify the TricValve® System identification number with the product label on the box.
7.	Remove the TricValve® System from the blister tray and inspect for any gross damage. Wipe the length of the catheter with moist (saline) gauze.
8.	Check that the TricValve® loaded capsule is slightly open to maximize flushing of the delivery system.  Figure 4: TricValve® capsule slightly open  Figure 5: Turn micro knob clockwise to unsheath device / open the capsule. If the TricValve® capsule is closed, use the micro knob on the handle of the delivery system and turn it clockwise (away from the operator) to open it until the tip is distal enough from the TricValve® capsule to allow for adequate flushing using sterile physiological saline solution.

9.	Fill tray with saline while submerging TricValve® Transcatheter Bicaval Valve System.
10.	Connect a 3-way stopcock to the tube flush port. Attach a 20 mL syringe filled with sterile physiological saline solution to the stopcock on the flush port. Flush the system with saline solution to ensure that the external tube is fully filled without bubbles. Repeat the step as needed with a minimum of 60 mL of saline solution and until all air bubbles are removed from the external tube. Keep the external tube filled with sterile physiological saline solution for a duration of minimum 2 minutes to keep TricValve® bioprosthesis hydrated (wet). It is crucial to follow these steps to ensure the appropriate hydration of the valve.  Figure 6: Attach syringe to stopcock and flush the external tube.
11.	Verify that no catheter leakage is observed during any of the flushing steps from the delivery system. If leakage is observed from any joints or section of the delivery system, replace, and use a new system.
12.	Close the TricValve® capsule using the micro knob and rotate it counterclockwise (towards the operator) on the delivery system handle until the micro knob can no longer be turned.  Figure 7: Turn micro knob counterclockwise to close the TricValve® capsule.  Figure 8: TricValve® capsule completely closed. Inspect that the TricValve® capsule is completely closed and there is no opening at top. If necessary, flush the delivery system again to remove any residual air bubbles by pointing the delivery system tip upwards.

13.	Attach the luer lock syringe filled with sterile physiological saline solution to the guidewire port at the end of the delivery system handle. Flush the guidewire path with saline solution.  Figure 9: Attach syringe to guidewire port and flush the guidewire path.
14.	Conduct a final inspection of the TricValve® delivery system for any damage or issues. Check the TricValve® loaded inside the delivery system capsule to make sure the frame is free of creases or folds beyond the second node from the inflow end. Ensure check is performed circumferentially around the entire bioprosthesis. The TricValve® System is now ready to use with TricValve® hydrated and air removed.

10.3. VASCULAR ACCESS

WARNING: Do not use excessive force and/or manipulation during advancement and positioning of the guidewire, as it may result in perforation/dissection of arteries, veins, and/or other cardiac structures. This may also result in cardiac arrhythmias and conduction disturbances.

Step	Procedure
1.	Prepare for femoral venous access using conventional percutaneous puncture technique. This can be achieved as per hospital protocol (either percutaneous access using Seldinger technique or surgical cutdown). Three access sites will be made for the procedure.
2.	Prepare 1 st venous access in right common femoral vein using 6 Fr introducer sheath. This access shall be used for device deployment.
3.	Prepare 2 nd venous access in left common femoral vein using 6 Fr introducer sheath. This access shall be used to mark the crossing of right pulmonary artery (rPA) with superior vena cava (SVC).

4.	Prepare 3 rd arterial access in left femoral artery using a 6F introducer sheath. This should be used for continuous blood pressure monitoring.
5.	Administer anticoagulant according to physician / clinical judgement. If Heparin is administered as an anticoagulant, check the activated clotting time (ACT) after initial bolus of heparin and recheck every 30 minutes thereafter. Maintain ACT $\geq$ 250 seconds. NOTE: Anticoagulant may be administered at any time prior to this point but avoid delaying beyond this point.

10.4. IMPLANTATION OF SVC VALVE

WARNING: Avoid excessive movement of the delivery system while executing the procedure to protect vasculature or cardiac structures. Avoid excessive rotation of the delivery system to maintain functionality of the delivery system.

Step	Procedure
1.	Prepare TricValve® Transcatheter SVC Valve System as per device preparation steps mentioned in Section 10.2. The TricValve® should be properly hydrated using sterile physiological saline solutions for more than 2 minutes.
2.	Place a 5 Fr arterial catheter in the left femoral artery for continuous invasive blood pressure monitoring on a hemodynamic monitoring system. Position the catheter under fluoroscopy at the proper location.
3.	Position a pulmonary catheter through the left femoral vein in the right pulmonary artery (rPA) to mark the crossing of the rPA with the SVC.
4.	Place a 6 Fr pigtail (multipurpose) catheter through Introducer sheath in right femoral vein and obtain an angiogram of the superior vena cava (SVC) with special attention to the confluence and the distance between the innominate vein ostium and SVC, at the level where the right branch of the pulmonary artery crosses. This is recommended for planning the deployment. CAUTION: Careful handling of multipurpose catheters may be necessary to cross the venous valves at entrance to subclavian or jugular vein.
5.	Exchange the pigtail catheter for a straight 0.035-inch extra-stiff guidewire with a soft tip and place it from the femoral vein access up to the subclavian or internal jugular vein. Establish proper guidewire path under fluoroscopy.
6.	Make a small incision at the access site and dilate the Venous access site in right femoral vein, as needed using step-up dilators to 27.5 Fr.
7.	Insert the TricValve® delivery system with SVC valve over the guidewire.






8.	Advance the TricValve® delivery system through the femoral vein and inferior vena cava (IVC) into the right atrium (RA) and superior vena cava (SVC).
9.	All the steps for TricValve® deployment must be performed under fluoroscopy guidance.
10.	Place the upper part of the SVC Valve in the confluence, with the belly of the SVC valve positioned above the right pulmonary artery (rPA) crossing. Note: Delivery system and guidewire position may be adjusted during valve placement to optimize implant position in the SVC.
11.	Using fluoroscopic guidance and echocardiography (optional), verify and confirm the delivery system position in the SVC.
12.	Under fluoroscopy guidance, the SVC valve can be recaptured and repositioned up to 50% deployment.  Figure 10: Location up to which SVC Valve may be recaptured. CAUTION: Once the valve capsule is retracted more than 50% to expose the TricValve®, the valve is unable to be retrieved or recaptured into the delivery system.
13.	The deployment of SVC valve shall be performed by slowly unsheathing the TricValve® device from the valve capsule at the correct position by turning the micro knob in clockwise direction (away from the operator) as shown in figure below:  Figure 11: Turn micro knob clockwise to unsheath device / open the capsule.
14.	If it is required to recapture the partially released TricValve®, turn the micro knob in counterclockwise direction until the valve capsule is completely closed as shown in figure below:  Figure 12: Partially released device from the delivery system



	 Figure 13: Rotate micro knob counterclockwise to recapture Device
15.	The SVC valve is a self-expanding device and shall take the shape of the anatomy on deployment. The sizing of the device ensures that the self-expanding valve fits properly in the superior vena cava.
16.	After satisfactory position is confirmed, slowly release the SVC valve from the delivery system using the micro knob and by turning it till it no longer turns. Under fluoroscopy observe the gradual release of the SVC valve in desired position. Double check that micro knobs no longer turn clockwise, and complete release is achieved.
17.	Gently rotate the delivery system in clockwise direction (<180°) and then back in counterclockwise direction (<180°) to confirm TricValve® device release from all the three anchors of delivery system.
18.	On confirmation of full deployment of TricValve® device, slowly withdraw the delivery system with open valve capsule till the tip has passed completely through the valve.
19.	Verify the delivery system – tip has passed through the valve using Fluoroscopy.
20.	Unlock the micro knob and retrieve the tip.  Figure 14: Unlock the micro knob and pull back of back handle of delivery system. Using your left hand, hold the red macro slider and slowly pull back the back-handle of the delivery system until the tip is fully




	retracted. The red macro slider must be kept in a fixed position during this process as shown below in Figure.  Figure 15: Tip retraction while keeping macro knob in fixed position.
21.	Close the TricValve® capsule and lock the micro knob as shown in figure below.  Figure 16: Lock the micro knob and close the TricValve® capsule
22.	Keeping the guidewire in position, slowly withdraw and remove the TricValve® delivery system, ensuring that the tip is safely retracted. Do not remove the guidewire and maintain its position in patient. CAUTION: Ensure the capsule is completely closed before catheter removal. If increased resistance is encountered when removing the catheter through the introducer sheath, do not force passage. Increased resistance may indicate a problem and forced passage may result in damage to the device and/or harm to the patient. If the cause of resistance cannot be determined or corrected, remove the catheter and introducer sheath as a single unit over the guidewire, and inspect the catheter and confirm that it is complete. CAUTION: Take care to safely withdraw the tip from valve and close the delivery system capsule, as failure to do so may impact valve functionality or lead to dislodgement of the valve. CAUTION: Failure to follow prescribed release steps may result in difficulty or inability to release implant, requiring additional intervention.
23.	Using fluoroscopy and echocardiography (optional) imaging, assess the final position and functionality of the SVC Valve.
24.	Perform pressure measurements to ensure correct function of SVC bioprosthesis. The next step shall be preparing and implanting IVC bioprosthesis.

10.5. IMPLANTATION OF IVC VALVE

WARNING: Avoid excessive movement of the delivery system while executing the procedure to protect vasculature or cardiac structures. Avoid excessive rotation of the delivery system to maintain functionality of the delivery system.

Step	Procedure
1.	Prepare TricValve® Transcatheter IVC Valve System as per device preparation steps mentioned in Section 10.2. The TricValve® should be properly hydrated using sterile physiological saline solutions for more than 2 minutes.
2.	Insert the TricValve® delivery system with IVC valve over the guidewire. CAUTION: All steps of implantation shall be performed under fluoroscopy imaging.
3.	All the steps for TricValve® deployment must be performed under fluoroscopy guidance.
4.	Advance the TricValve® delivery system through the femoral vein into the inferior vena cava (IVC). Position the IVC Valve at the height of diaphragm with the skirt visible just above the hepatic vein inflow.
5.	Use of a 5 Fr pigtail (multipurpose) catheter placed in upper hepatic vein may be used for landmarking purposes.
6.	Align the constrained segment of the IVC Valve stent frame with the Cavo-atrial junction by carefully pulling back the TricValve® delivery system. Note: Delivery system and guidewire position may be adjusted during valve placement to optimize implant position in the IVC.
7.	Using fluoroscopic guidance and echocardiography (optional), verify the delivery system position alignment with Cavo-atrial junction in the IVC. Removal of 5 Fr pigtail catheter placed in upper hepatic vein can take place up to 50% of IVC valve deployment. CAUTION: With a safety margin of 5 mm, take precautions for a low or high valve position in IVC, as it can cause either hepatic vein obstruction or a paravalvular leak.
8.	The deployment of IVC valve shall be performed by slowly unsheathing the TricValve® device from the valve capsule at the correct position by turning the micro knob in clockwise direction (away from the operator), as show in Figure below.  Figure 17: Turn micro knob clockwise to unsheath device / open the capsule.

9.	Under fluoroscopy guidance, the IVC valve can be recaptured and repositioned up to 50% deployment as shown in figure below.  Figure 18: Location up to which IVC Valve may be recaptured. CAUTION: Once the valve capsule is retracted more than 50% to expose the TricValve®, the valve is unable to be retrieved or recaptured into the delivery system.
10.	If it is required to recapture the partially released TricValve® and reposition it, turn the micro knob in counterclockwise direction until the valve is captured and valve capsule is completely closed as shown in figure below.  Figure 19: Partially released device from the delivery system  Figure 20: Rotate micro knob counterclockwise to recapture device
11.	The IVC valve is a self-expanding device and shall take the shape of the anatomy on deployment. The correct sizing of the device ensures that the self-expanding valve fits properly in the inferior vena cava.
12.	After satisfactory position is confirmed, slowly release the IVC valve from the delivery system using the micro knob and by turning it till it no longer turns. Under fluoroscopy observe the gradual release of the IVC valve in desired position. Double check that micro knob, no longer turns clockwise, and complete release is achieved.

13.	Gently rotate the delivery system in clockwise direction (<180°) and then back in counterclockwise direction (<180°) to confirm TricValve® device release from all the three anchors of delivery system.
14.	On confirmation of full deployment of TricValve® device, slowly withdraw the delivery system with open valve capsule till the tip has passed completely through the valve.
15.	Unlock the micro knob and retrieve the tip. Verify the delivery system – tip has passed through the valve using fluoroscopy.  Figure 21: Unlock the micro knob and pull back of back handle of delivery system. Using your left hand, hold the red macro slider and slowly pull back the back-handle of the delivery system until the tip is fully retracted. The red macro slider must be kept in a fixed position during this process as shown below in Figure.  Figure 2: Tip retraction while keeping macro knob in fixed position.
16.	Close the TricValve® capsule and lock the micro knob as shown in figure below.  Figure 3: Lock the micro knob and close the TricValve® capsule
17.	Keeping the guidewire in position, slowly withdraw and remove the TricValve® delivery system, ensuring that the tip is safely retracted. Do not remove the guidewire and maintain its position in patient. CAUTION: Ensure the capsule is completely closed before catheter removal. If increased resistance is encountered when removing the catheter through the introducer sheath, do not force passage. Increased resistance may indicate a problem and forced passage may result in damage to the device and/or harm to the patient. If the cause of resistance cannot be determined or corrected, remove the catheter and introducer sheath as a single unit over

	the guidewire, and inspect the catheter and confirm that it is complete. CAUTION: Take care to safely withdraw the tip from valve and close the delivery system capsule, as failure to do so may impact valve functionality or lead to dislodgement of the valve. CAUTION: Failure to follow prescribed release steps may result in difficulty or inability to release implant, requiring additional intervention.
18.	Using fluoroscopy and echocardiography (optional) imaging, assess the final position and functionality of the IVC Valve.
19.	Perform pressure measurements to ensure correct function of IVC bioprosthesis.

10.6. POST DEPLOYMENT


Step	Procedure
1.	Place a 6 Fr pigtail (multipurpose) catheter through the right femoral vein on guidewire and obtain an angiogram of the implants. This is not mandatory but highly recommended for visualizing the correct functioning of the TricValve® bioprosthesis in Superior and Inferior vena cava.
2.	Remove the catheter inserted and Introducer sheath from left and right access sites as per standard hospital protocol.
3.	Perform femoral access site closure as applicable using standard percutaneous intervention techniques. CAUTION: Use of Percutaneous closure devices are recommended whenever sheath size of 6 Fr or larger is used in a percutaneous interventional procedure. This is to reduce risk of hematoma and aid in the hemostasis of vascular access site.
4.	The final position of TricValve® bioprostheses in Superior and Inferior vena cava shall resemble the final implantation result image as shown below: 

	Figure 4: Final implant position in Superior and Inferior vena cava.
5.	Perform Transesophageal Echocardiography (TEE) / Transthoracic Echocardiography (TTE) immediately after the TricValve® implantation procedure, if needed.
6.	Administer anti-coagulant and / or anti-platelet therapy, as required according to hospital protocol.

11 HOW SUPPLIED

11.1. STERILIZATION AND PACKAGING

The TricValve® Transcatheter Bicaval Valve system is provided sterile and non-pyrogenic in individual device packaging (two separate boxes, one for SVC and another one for IVC device). The bioprosthesis is already loaded inside each delivery system. The TricValve® Transcatheter Bicaval Valve system is packaged in a cardboard box containing double Tyvek® pouch with the device. The TricValve® Transcatheter Bicaval Valve System is terminally sterilized using ethylene oxide (EO) gas.

WARNING: The outer surface of the sealed 1st Tyvek® pouch is not sterile and must not be placed in the sterile field.

11.2. STORAGE

The TricValve® Transcatheter Bicaval Valve system should be stored in a clean dry area that is contamination free and at temperature between 8°C to 25°C [46.4°F to 77°F].

CAUTION: Avoid exposing the Device package to extreme temperature fluctuations and direct contact with sunlight. Avoid freezing temperatures as it can damage the Bioprostheses. Storage at elevated temperatures may potentially damage the polymeric components and adhesives thus compromising the product performance.

12 MAGNETIC RESONANCE IMAGING (MRI) SAFETY INFORMATION

MR Conditional

Non-clinical testing has demonstrated that the TricValve® Transcatheter Bicaval Valve is MR Conditional. A patient with this device can be safely scanned in an MR system meeting the following conditions, failure to follow these may result in injury.

- Static magnetic field of 1.5 T or 3.0 T (Tesla).
- Maximum spatial gradient magnetic field of ≤ 4000 gauss/cm (40 T/m) or less.
- Maximum MR system-reported, whole body averages specific absorption rate (SAR) of 2.0 W/kg
- Normal operating mode of the MR System for both gradients and SAR

Under the scan conditions defined above, the TricValve® implant is expected to produce a maximum temperature rise of less than 4 °C after 15 minutes of continuous scanning.

In non-clinical testing, the image artifact caused by the TricValve® implant extends approximately 1.30 cm from the device in Spin Echo (SE) and 1.02 cm in Gradient Echo (GE) sequence, when imaged with 3.0T MRI system.

NOTE: The TricValve® Transcatheter Bicaval Valve System has not been evaluated in MR systems other than 1.5 or 3.0 Tesla.

13 EXPECTED LIFETIME OF DEVICE

The TricValve® Transcatheter Bicaval Valve System has been subjected to rigorous pre-clinical durability testing per the valve testing requirements (ISO 5840 standard). The valves were successfully tested to 5 years of simulated wear. Actual lifetime performance in humans depends on many factors and varies from patient to patient. Specific activities or conditions that could shorten or lengthen the device's lifetime have not been established.

14 PATIENT INFORMATION

A Patient implant card (IC) is provided with each TricValve® Transcatheter Bicaval Valve System package. After implantation, please complete all requested information and provide the implant card to the patient. The serial number is found on the TricValve® package. This implant card allows patients to inform healthcare providers what type of implant they have when they seek care.

15 RECOVERED IMPLANTS AND DEVICE DISPOSAL

P+F Products + Features GmbH is interested in obtaining explanted TricValve® Transcatheter Bicaval Valve and its delivery system for analysis. A written report summarizing our findings will be provided upon completion of evaluation. Please contact P+F Products+Features GmbH for return of the recovered valves.

If you decide to return any of the devices, please follow the following instructions:

- Unopened package with sterile barrier intact: If the pouches have not been opened, return the device in its original packaging.
- Package opened but not implanted: If a pouch is opened, the device is no longer sterile. Please return the device with its original packaging.
- Explanted implant: The explanted implant should be placed into a suitable histological fixative such as 10% formalin or 2% glutaraldehyde. Refrigeration is not necessary under these circumstances. Contact company personnel to request return to P+F Products+Features GmbH.

15.1. DISPOSAL OF DELIVERY SYSTEM

Use universal precautions for biohazards and sharps to avoid user injury. Used single use devices (including all devices that come in direct contact with patient) should be handled and disposed of in accordance with institutional

guidelines on biohazardous materials and hospital waste to avoid possible cross-contamination.

The Used delivery system may be disposed in the same manner as the hospital waste and biohazardous materials are handled. There are no special or unusual risks related to the disposal of delivery system which are required to be followed.

16 SUMMARY OF CLINICAL EXPERIENCE

A high-level summary of all clinical experience is provided below. A detailed summary of all clinical experience can be found in the SSCP. The SSCP (Summary of Safety and Clinical Performance) has been adapted in accordance with the clinical evaluation assessment by the Notified Body on which the CE certification has been granted. The SSCP contains a relevant summary of the same information.

Conformity of the TricValve® Transcatheter Bicaval Valve System of the General Safety and Performance Requirements (GSPR) for safety (MDR GSPR 1), performance (MDR GSPR 1), acceptability of side-effects (MDR GSPR 8), usability (MDR GSPR 5), device lifetime (MDR GSPR 6), acceptable benefit-risk profile (MDR GSPR 8) has been established for the labelled indications. After the launch of the European Database on Medical Devices/Eudamed, refer to <https://ec.europa.eu/tools/eudamed> for a SSCP for this medical device.

16.1. TRICUS STUDY

The TRICUS Study is a prospective, single-arm, single centre study, planned as an early feasibility / first-in-human study, involving 9 patients with severe symptomatic TR at extremely high operative risk as assessed by Heart Team. The Study showed a high rate of device, procedural and clinical success. The implantation success of 89% was observed with 6-month survival rate of 100% (9 out of 9 patients alive). At 6-months no MAE (major adverse event) were reported. Immediate hemodynamic success was measured by reduction of the venous pressure IVC: v-wave (mmHg) from 21.9 ± 6.4 mmHg before IMD implantation to 16.6 ± 4.7 mmHg after the procedure. The improvement of clinical parameters showed an improvement in quality of life, exercise capacity, and functional status. The Quality of life measured by the mean KCCQ (Kansas City Cardiomyopathy Questionnaire) showed improvement at 6 months. The Exercise capacity of patients as assessed by 6-minute walk test (6 MWT) improved at 6 months.

16.2. TRICUS EURO STUDY

The TRICUS Euro pivotal study is a prospective, multi-national, multi-centre, non-blinded, non-randomized single-arm clinical trial conducted to evaluate the safety and effectiveness of TricValve® Transcatheter Bicaval Valve System in the Superior and Inferior vena cava in patients with severe Tricuspid Regurgitation (TR). The study represents experience of TricValve® System in NYHA (New York Heart Association) functional class III or IV severe TR patients, optimally medicated and ineligible for open heart surgery, with significant caval backflow. The Clinical Study enrolled a total of 35 patients from 12 institutions in Spain and Austria. All enrolled study patients shall be followed

at 30 days, 3 months, 6 months, 1 year, and annually for 5 years. The implantation (procedure) success was 91% with a 6-month survival rate of 91% (32 out of 35 patients alive). At 6 months, TricValve® System treatment reduced heart failure symptoms in 76% of the patients (improvement in NYHA class by at least 1 grade); quality of life (KCCQ score) increased from 41 ± 22 points to 56 ± 24 points.

16.3. CLINICAL DATA ANALYSIS

Data below represents the results of TRICUS and TRICUS EURO Clinical study outcomes at 1 year for all the patients enrolled in the studies. Both the studies utilized an independent CT core lab and Echocardiography core lab, to assess patient data. An independent CEC adjudicated safety events and a DSMB independently reviewed aggregate safety data and evaluated trends of adverse events and their effect on trial conduct and device risk assessment. A total of 44 study participants were evaluated.

	6 mo	1 y	Overall 1 y (%)
Major adverse events			
All-cause death	3	0	6.8
Myocardial infarction	0	0	0
Tricuspid valve surgery	1	0	2.3
Cardiac tamponade	0	0	0
Major bleeding ^a	7	2	20.4
Stroke	3	1	9
Serious adverse events			
Heart failure rehospitalization	9	4	29.5
Right heart thrombi	4	2	13.6
Paravalvular leak	3	0	6.8

^aMajor bleeding Valve Academic Research Consortium-2 criteria: overt bleeding either associated with a drop in the hemoglobin level of at least 3 g/dL or requiring transfusion of 2 to 3 units of whole blood/red blood cells, or causing hospitalization or permanent injury, or requiring surgery and does not meet criteria of life-threatening or disabling bleeding.

Figure 5: Major and serious adverse events at 1-year follow-up.

Key Demographics and Baseline Characteristics:
The Enrolled patients had a mean age of 76.2 ± 7.5 years with 82% female patients. All patients were highly symptomatic (86.4% in NYHA functional class III and 13.6% in NYHA functional class IV). Most of the patients had a history of prior valve interventions (72%) and significant comorbidities like Atrial fibrillation (93%), prior pacemaker leads (20.5%) and renal dysfunction (65.9%). The overall EUROSCORE II was 5.6% ± 3.6% and TRISCORE was 5.4 ± 1.2 for the group. All patients presented with symptoms of right heart failure, such as lower limb Edema, ascites, or pleural effusion, and the most common right heart failure signs were fatigue, shortness of breath, and exercise intolerance. Hepatic vein backflow was present in 86.5% of the sample. Mean systolic pulmonary artery pressure estimated by echocardiography was 43 ± 9 mm Hg.

Safety. Clinical Studies met the primary safety endpoint. The major adverse event (MAE) rate at overall 1 year for TricValve® patients was 29.5%, which is less than the performance level of 70%, based on surgical predicates. The CEC-adjudicated MAEs at 1 year are shown below.


Figure 6: KCCQ-12 Score.


Figure 7: NYHA Functional Class.

There were 3 deaths during the 6-month follow-up (1 due to cardiovascular cause) with no further deaths at 1 year follow-up. Following implantation of TricValve®, there was only one incidence of SVC valve embolization which occurred day after the implantation and required emergent cardiac surgery. No cases of myocardial infarction nor cardiac tamponade were detected. No evident harm to patients have been identified for patients with mild valvular and paravalvular leakage at 1 year.

Primary Effectiveness Endpoints: The primary effectiveness endpoint was clinical improvement evaluated by the composite of: change in quality of life (QOL) measured with the 12-item Kansas City Cardiomyopathy Questionnaire (KCCQ-12), with a large improvement defined as an increase of ≥15 points from baseline to 1-year follow-up; improvement in NYHA functional class to I or II at 1-year follow-up; or change in functional exercise capacity measured with the 6MWT distance (≥ 40m).


Figure 8: 6-Minute Walk Test.

The study demonstrated that TricValve® Bicaval Valve System is more likely to improve symptomatic and functional outcomes. The present findings confirm these beneficial effects remain stable and significant out to 1-year in the combines TRICUS cohorts.

Secondary Endpoints: The secondary endpoints were primarily safety related and are as follows: 1) freedom from major adverse events including death, acute myocardial infarction, tricuspid valve surgery, cardiac tamponade, stroke, or major bleeding according to Valve Academic Research Consortium-3 (VARC-3) criteria at 1-year follow-up; 2) freedom from heart failure rehospitalizations or serious adverse events related to the device at 1-year follow-up; 3) changes in right heart dimensions as measured by assorted echocardiographic parameters at baseline and 3-month, 6-month, and 1-year follow-up) improvement in systemic venous congestion measured by echocardiographic and laboratory test parameters (hepatic vein backflow and N-terminal pro-B-type natriuretic peptide); and 5) changes in renal and hepatic function measured by baseline and 1-year follow-up laboratory tests (creatinine, glomerular filtration rate, alanine aminotransferase, and aspartate aminotransferase).

Results: Primary combined effectiveness endpoint of clinical improvement was achieved in 95.5% (42) patients. Improvement in QOL as measured by an increase of at least 15 points in the KCCQ-12 score was seen in 56.4% of the patients (Figure 26). A significant improvement in NYHA functional class was observed in 62.2% of patients in NYHA functional class I or II at 1-year follow-up. Additionally, in Exercise capacity evaluated and showed large improvement in 40% of patients in 6MWT distance ($\geq 40m$) at 1 year after implantation (Figure: 28).

Results: Major adverse events at 6-month and 1-year follow-up are shown above in Figure 25. There were 3 deaths during the 6-month follow-up (1 due to cardiovascular cause), with no further deaths at 1-year follow-up. No cases of myocardial infarction nor cardiac tamponade were detected, but 1 patient in the TRICUS first-in-human cohort required emergent cardiac surgery the day after

the implantation due to embolization of the SVC prosthesis during deployment. The Kaplan-Meier survival curve is presented in Figure 29(a).

Major bleedings were essentially gastrointestinal hemorrhages (9.0% of the population) and puncture site complications (6.8%), from which only 1 patient needed surgical vascular repair. Additionally, 4 (9.0%) patients experienced a stroke during follow-up. The heart failure rehospitalization rate was 29.5% at 1-year follow-up and occurred more often in those patients with prior mitral percutaneous repair but was not related to prior mitral surgery. For cases of right heart chamber thrombi in 6 cases, it was incidental findings in scheduled follow-up and was medically managed with no clinical consequences. Freedom from re-Hospitalization Kaplan-Meier curve is presented in Figure 29(b).


Figure 9: Kaplan - Meier Curve (a) 1-Year survival (b) Freedom from Hospitalization.

One-year echocardiographic parameters showed no significant differences on right ventricular function as assessed by TAPSE were detected, but there was a significant reduction in tricuspid annular diameter at 1-year follow-up. Hepatic vein backflow was abolished in 63.8% of patients at hospital discharge post-CAVI and remained unchanged throughout the rest of the follow-up period, resulting in reduced congestion symptoms and less reliance on diuretic therapy. Regarding other laboratory test parameters, no significant changes were found in any of the studied variables.

The vast majority (95.0%) of TricValve recipients demonstrated some objective evidence of significant functional improvement, with 56.0% demonstrating a significantly large improvement in the KCCQ-12 score of ≥ 15 points from baseline. The heart failure hospitalization rate at 1 year was 29.5% with an all-cause mortality rate of 6.8%. The immediate abolition of reverse hepatic vein flow on Doppler imaging post-CAVI significantly associated with a large (delta KCCQ-12 score of ≥ 15 points from baseline) functional improvement at 1 year.

16.4. ADDITIONAL STUDIES:

TRICUS REGISTRY - Retrospective Prospective Multicentric Clinical Follow up of Patients With Severe Tricuspid Regurgitation After Being Treated With the TricValve® Transcatheter Bicaval Valves System. Post market clinical follow-up study is currently ongoing and has yet to reach their primary endpoints; therefore, results are not presented here.

16.5. CLINICAL CONCLUSION:

The Data collected in the TRICUS and TRICUS EURO Clinical trials demonstrates that TricValve® Transcatheter Bicaval Valve System is a safe and effective treatment option for patient with severe symptomatic Tricuspid Regurgitation / Right heart failure population who currently have limited to no treatment options. The benefits of TricValve® System outweigh the risks of its use and provide a reasonable assurance of safety and effectiveness. The Clinical studies of TricValve® Transcatheter Bicaval Valve System met the primary safety and effectiveness endpoints, demonstrating (1) an acceptable safety profile for TricValve® Transcatheter Bicaval Valve System, (2) effective improvement in KCCQ-12 (quality of life), NYHA function classification and 6MWT (exercise capacity) (3) reduction or abolition of reverse hepatic backflow, and (4) increase in survival rates with decreased rehospitalization rates.

Long term follow-up data will be collected as part of 5-year follow-up for patients enrolled in both clinical trials to confirm safety and effectiveness of TricValve® Transcatheter Bicaval Valve System.

17 CLINICAL BENEFITS

The clinical benefits of the TricValve® Transcatheter Bicaval Valve System for the treatment of severe tricuspid regurgitation (TR) include the following:

- Offers a safe treatment option for patients with severe tricuspid regurgitation.
- Enables a minimally invasive tricuspid valve regurgitation treatment.
- Simple procedure of implantation in Vena Cava without disturbing native tricuspid valve.
- Enables clinical improvements in quality of life, cardiac function, and exercise capacity.
- Mechanical ventilation and transesophageal echocardiography not mandatory required.
- High technical success in implantation of device in Vena Cava.
- Significant reduction in Hepatic Vein congestion.
- Independent of any pre-existing native tricuspid valve deformation.
- Independent of any pre-existing pacemaker leads or prior valve surgery.

- Not disturbing the native tricuspid valve, allowing further treatment options for the future.

Treatment with the TricValve® Transcatheter Bicaval Valve System may lead to elimination or a significant reduction of TR, which can result in clinically meaningful improvements in quality of life, cardiac function, and exercise capacity. The clinical improvement depends on the seriousness and reversibility of the tricuspid regurgitation-associated disease.

18 BASIC UDI-DI (Unique Device Identification – Device Identifier)

The Basic UDI-DI is the access key for device-related information entered in the Eudamed. The Basic UDI-DI for the TricValve® Transcatheter Bicaval Valve System can be used to locate the SSCP. The following table 3 contains the Basic UDI-DIs for the TricValve® Transcatheter Bicaval Valve System:

Table 3 Unique Device Identification of TricValve® variants

Medical Device Name	Basic UDI-DI
TricValve® Transcatheter Bicaval Valve System-SVC 25	04260642340001
TricValve® Transcatheter Bicaval Valve System-SVC 29	04260642340018
TricValve® Transcatheter Bicaval Valve System-IVC 31	04260642340025
TricValve® Transcatheter Bicaval Valve System-IVC 35	04260642340032

Disclaimer (Use & Responsibility)

The TricValve® Bicaval Transcatheter Valve System (Device) is a Class III implantable medical device intended for use exclusively by qualified and appropriately trained healthcare professionals in accordance with the instructions, indications, contraindications, warnings, precautions, and training requirements described in this Instructions for Use (IFU). This IFU is intended to support the safe and effective use of the Device and does not replace the clinical judgment of the treating physician. Patient selection, procedural planning, implantation technique, and post-implantation management remain the sole responsibility of the healthcare professional. The Device must be used strictly in accordance with its intended purpose and the instructions provided by P+F Products + Features GmbH. Use outside the scope of the intended purpose, including off-label use, reuse, unauthorized modification, improper handling or storage, or use by unqualified personnel, may result in serious injury or death and is not covered by the manufacturer. As with all implantable heart valve systems, use of the Device involves inherent risks, including the risk of serious injury or death. Failure to follow the instructions and warnings contained in this IFU may compromise device performance and patient safety. To the extent permitted by applicable law, P+F Products + Features GmbH shall not be responsible for adverse events arising from improper use of the Device, deviation from this IFU, or failure to comply with applicable clinical standards of care. The application of this Disclaimer is subject to applicable local laws and regulations in the country of use. Nothing in this section is intended to limit or exclude any responsibility or liability that cannot be excluded under mandatory applicable law, nor does it affect the manufacturer's obligations relating to device safety, vigilance, or corrective actions under applicable medical device regulations.

Limited Warranty

P+F Products + Features GmbH warrants that the TricValve® Bicaval Transcatheter Valve System has been designed, manufactured, and tested in accordance with applicable regulatory requirements and quality standards and

is free from defects in materials and workmanship at the time the Device is released by the manufacturer. Due to the nature of implantable medical devices and patient-specific anatomical and clinical factors, no guarantee can be given that the Device will perform satisfactorily in all patients or clinical situations. This limited warranty applies only when the Device is used in accordance with its intended purpose and this IFU. The warranty does not apply in cases of misuse, off-label use, improper storage or handling, unauthorized modification, reuse, or failure to comply with the instructions and warnings provided. To the extent permitted by applicable law, no additional express or implied warranties are provided beyond those expressly stated herein. Any remedies under this limited warranty are limited to those required under applicable law. This limited warranty does not limit or exclude any rights of patients or users under mandatory local, national, or international law, including applicable product liability and patient protection legislation. The application of this Limited Warranty is subject to applicable local laws and regulations in the country of use. If any provision of this Limited Warranty is held to be invalid or unenforceable by a competent authority, the validity of the remaining provisions shall not be affected.

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Box plot showing GMWT (m) for five groups. The y-axis represents GMWT (m) from 0 to 500. The x-axis categories are: Pradinis rodiklis (n=44), 30 dienų (n=30), 3 mėnesiai (n=29), 6 mėnesiai (n=34), and 1 metai (n=32). The plot shows a general upward trend in median GMWT over time. Individual data points are overlaid on the box plots. Statistical significance is indicated by p-values: P=0.902 for the first two groups, P=0.657 for the first three, and P=0.290 for the first four.

Group	n	Median GMWT (m)	Q1 (m)	Q3 (m)	Min (m)	Max (m)
Pradinis rodiklis	44	253.4	160.0	302.5	~60	~400
30 dienų	30	235.5	168.0	324.0	~80	~460
3 mėnesiai	29	267.0	160.0	309.0	~10	~460
6 mėnesiai	34	275.0	200.0	319.0	~60	~460
1 metai	32	284.0	197.5	341.0	~80	~460

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INFORMATION ON MATERIAL AND SUBSTANCES INCORPORATED

For the TricValve® Transcatheter Bicaval Valve System, the following table 4 shows the information on the materials and substances incorporated in the Medical Device – TricValve® Valve:

Table: Qualitative and Quantitative information of material:

Substance	CAS Number	Model Mass Range (mg)
Nickel	7440-02-0	1113 – 1504
Titanium	7440-32-6	835 – 1128
Cobalt	7440-48-4	0.00 – 1.320
Iron	7439-89-6	0.00 – 1.320
Niobium	7440-03-1	0.00 – 0.660
Chromium	7440-47-3	0.00 – 0.262
Copper	7440-50-8	0.00 – 0.262
Carbon	7440-44-0	0.00 – 1.060
Hydrogen	1333-74-0	0.00 – 0.131
Nitrogen	7727-37-9	0.00 – 0.131
Oxygen	7782-44-7	0.00 – 1.060
Collagens, Bovine Pericardium	NA	180 – 570
Polytetrafluoroethylene (PTFE)	9002-84-0	180 – 220
Polyester Fibers (PET fibers)	80595-68-2	405 – 570
Propane-1,2,3-triol	56-81-5	0.00 – 0.008
Androstan-17-one,3-ethyl-3-hydroxy-, (5.alpha.)-	57344-99-7	0.0341 – 0.0683
Carbonic acid, methyl ester,[(E,E)-3,7,11-trimethyl- 2,6,10-dodecatrien-1-yl] ester	1000164-38-7	0.0088 – 0.0177
26-Nor-5-cholesten-3.beta.-ol-25-one	7494-34-0	0.0219 – 0.0439
Supraene	7683-64-9	0.0048 – 0.0096
Cholesterol	57-88-5	0.0866 – 0.0958
1,4-p-Menthadien-7-al	22580-90-1	0.0005 – 0.0011
1-Isopropyl-2,2- dimethylpropylideneamine	1000210-05-3	0.0006 – 0.0011
Heptane, 2,4-dimethyl-	2213-23-2	0.0008 – 0.0016
Irganox PS800	123-28-4	0.0005 – 0.0011
Irganox 1010	6683-19-8	0.0049 – 0.0099
Oleamide	301-02-0	0.0004 – 0.0008
Erucamide - (Z)-Docos-13-enamide	112-84-5	0.0008 – 0.0017

P+F Products + Features GmbH
Landsberger Strasse 155, Building 3, 80687 Munich, Germany
Telephone: +49 89 255 448 400
Email : info@productsandfeatures.com
www.productsandfeatures.com

11-2025
IFU-HV01-02. V0
Copyright 2025, P+F Products+Features GmbH